


Inputs and Outputs

LabView			
Inputs		Outputs	
Neural Decode	MATLAB	Decode Method	MATLAB
Trained Transfer Matrix/Information	MATLAB	Training Task/Trial Information	MATLAB
Residual Limb Movements	Vicon	Task/Trial Information	Unity
Hand Movements	Unity	Vicon Control	Vicon
Audio	Microphone	Data Streams/Markers	Hard Drive

MATLAB			
Inputs		Outputs	
Training Task/Trial Information	LabView	Trained Transfer Matrix/Information	LabView
Decode Method	LabView	Neural Decode	LabView
Hand Movements	Unity	Neural Decode	Unity
Streaming Neural Data	TDT PO8e Card		

Unity			
Inputs		Outputs	
Hand/Finger Movements	Leap Motion	Hand/Finger Movements	MATLAB
Residual Limb Movements	Vicon	Hand/Finger Movements	LabView
Head Movements	Oculus Rift	Head Movements	Hard Drive
Neural Decode	MATLAB	Oculus Rift Video	Hard Drive
Task/Trial Information	LabView		

LabView Experiment Outline



- Initialize and start accessory programs
 - Via PowerShell, start MATLAB, Unity, etc...
 - Ping software to ensure proper connections.
- Equip the patient with the Oculus Rift, Vicon stickers, and patient cable
- Ask the patient to move their arms until Leap Identifies their intact hand
- Ask patient to move to a rest position for 10 sec (record this base or rest position for 'Return to Rest' cue)
 - Ask patient to mimic phantom perceived positions with intact hand for 10 sec (record this base or rest position for 'Return to Rest' cue)
- Start Tasks
 - Which finger(s) or hand synergies
 - Which finger task (spherical or torus)
 - Save timestamps
 - Update Unity and MATLAB
- Trial Params/Information:**
 - Task Duration
 - Inter-task interval
- Sphere Task**
 - Finger(s) or synergies
 - Hold Time
 - Success Radius
 - Return to Rest cue
- Torus Task**
 - Finger
 - Success Radius
 - Rotation Number
 - Target (blue sphere)
 - Target Rotation Direction
 - Target Speed (Ramp on/Ramp off)
 - Cursor Position (Fingertip center)
 - Torus Color (Red or green based on the success radius)
 - Return to Rest cue
- Save Streams
 - TimeStamps
 - Markers (success, failure, trial start, trial information/parameters)
 - Data Streams (Neural Decode, Finger Movements)

Task/Trial Summary**Patient Rest**

- Ask patient to be in a 'relaxed' or base position (Leap must be able to see all 5 fingers)
- 10 seconds hold time

Mimic Phantom Hand

- Ask patient to mimic the perceived position of their phantom hand with their intact hand
- 10 seconds hold time

Sphere Task

- Ask patient to move the indicated finger(s) or hand to the cued position
- Red colored spheres will appear for the selected finger(s) and the patient must move their fingers to that position.
 - Spheres will become green when the patient's fingertip(s) are within some success radius of the sphere center(s).
 - Patients will be cued to return to a rest or base position (positions gathered from the Patient Rest task).

Torus Task

- Ask patient to follow a blue target sphere around a torus
- A red torus will appear around the selected finger.
 - A blue target sphere will appear on the torus.
 - Patient will move their finger to the blue target sphere.
 - Target sphere will start moving (clockwise or counter-clockwise), and will slowly ramp up speed.
 - Torus will change from red to green when the patient's finger is within some success radius of the blue target sphere.
 - The blue target sphere will rotate around the torus one or more times, at a set speed (with a ramp up and a ramp down).
 - The blue target sphere will appear at the same position each trial. The blue target sphere will also stop at that position.


Arm Task?

- Ask patient to move arms from resting on the table to held in the air
- Perform the Sphere Task or Torus Task

Wrist Task?

- Ask patient to move wrist at different angles
- Perform the Sphere Task or Torus Task or Arm Task?

UDP Connection Scheme



Unity -> LabView

- Data: Hand positions
- IP: 192.198.0.20
- Port: 9001

LabView -> Unity

- Data: Target information
- IP: 192.198.0.10
- Port: 9000

MATLAB -> LabView

- Data: Neural Decode/Trained Matrix/Transfer Matrix/...
- IP: 127.0.0.1
- Port: 9011

LabView -> MATLAB

- Data: Task Information, Success information, (For training)
- IP: 127.0.0.1
- Port: 9010

Unity -> MATLAB

- Data: Hand Positions
- IP: 192.198.0.20
- Port: 9021

MATLAB -> Unity

- Data: Hand Positions
- IP: 192.198.0.10
- Port: 9020

Vicon -> LabView

- Data: Arm Positions
- IP: 192.198.0.20
- Port: 9031

Vicon -> Unity

- Data: Hand Positions
- IP: 192.198.0.20
- Port: 9030

LabView -> Vicon

- Data: ???
- IP: 192.198.0.30
- Port: 9032

Packet Construction

LabView -> Unity

Single Finger Trial Example

Trial Start:

'T:0;D:0;'

- Trial: 0 (Sphere); Display: 0 (off);

Trial Cue:

'F:0;Q:[s,i,j,k][s,i,j,k][s,i,j,k];C:R;'

- Finger: 0 (Thumb); Quaternion: [s,i,j,k][s,i,j,k][s,i,j,k]; Color: Red;

Display Trial:

'F:0;D:1;'

- Display: 1 (on);

Success Visual:

'F:0;C:G;'

- Finger: 0 (Thumb); Color: Green;

Return to Rest Cue:

'F:0;Q:[s,i,j,k][s,i,j,k][s,i,j,k];C:R;'

- Finger: 0 (Thumb); Quaternion: [s,i,j,k][s,i,j,k][s,i,j,k]; Color: Red;

Display Trial:

'D:0;'

- Display: 0 (off);

Multi Finger Trial Example

Trial Start:

'T:0;D:0;'

- Trial: 0 (Sphere); Display: 0 (off);

Trial Cue

'F:012;Q:[s,i,j,k][s,i,j,k][s,i,j,k][s,i,j,k][s,i,j,k][s,i,j,k][s,i,j,k][s,i,j,k];C:RRR;'

- Finger: 0,1,2 (Thumb, Index, Middle); Quaternion: [s,i,j,k][s,i,j,k][s,i,j,k] (First Finger), [s,i,j,k][s,i,j,k][s,i,j,k] (Second Finger), [s,i,j,k][s,i,j,k][s,i,j,k] (Third Finger); Color: Red, Red, Red;

Display Trial:

'D:1;'

- Display: 1 (on);

Partial Success Visual

'F:01;C:GG;'

'F:2;C:G;'

Return to Rest Cue

'F:012;Q:[s,i,j,k][s,i,j,k][s,i,j,k][s,i,j,k][s,i,j,k][s,i,j,k][s,i,j,k][s,i,j,k];C:RRR;'

Display Trial:

'D:0;'

- Display: 0 (off);

Torus Task Example

Trial Start:

'T:1;D:0;'

- Trial: 1 (Torus); Display: 0 (off);

Trial Cue:

'F:1;Q:[s,i,j,k][s,i,j,k][s,i,j,k];C:R;'

- Finger: 1 (Index); Quaternion: [s,i,j,k][s,i,j,k][s,i,j,k]; Color: Red;

Display Trial:

'D:1;'

- Display: 1 (on);

Success Visual:

'F:0;C:G;'

- Finger: 1 (Index); Color: Green;

Return to Rest Cue:

'F:0;Q:[s,i,j,k][s,i,j,k][s,i,j,k];C:R;'

- Finger: 1 (Index); Quaternion: [s,i,j,k][s,i,j,k][s,i,j,k]; Color: Red;

Display Trial:

'D:0;'

- Display: 0 (off);

Unity -> LabView/MATLAB

Hand/Finger Positions

'F:012345;Q:[s,i,j,k][s,i,j,k][s,i,j,k];'